

PREMIXES AND FEED ADDITIVES

SANI HORSE POTENCY



COMPOSITION - each one Kg contains:

ITEMS	AMOUNT	UNIT
Lipoic Acid (Tioctic Acid)	10 000	mg
Creatine	200 000	mg
Magnesium Sulphate Heptahydrate	3 000	mg
Additives (per Kg)		
Vitamins and provitamins		
Vitamin B1 (3a821 Thiamine Mononitrate)	12 000	
Vitamin B2 (Riboflavin-5 Phosphate)	1800	
Calcium D-Pantothenate (3a841)	800	
Vitamin B6/Pyridoxine Hydrochloride (3a831)	1000	
Vitamin B12 (Cyanocobalamin)	15	
Vitamins, provitamins and chemically well-defined sul	bstances having sir	nilar effect
Choline Chloride (3a890)	2 000	mg
Chemically well-defined substances having a similar bi	ological effect of vi	itamins
Niacinamide (3a315 Vitamin PP)	9 000	mg
L-Carnitine (3a910)	40 000	mg
Betaine HCl/Betaine Hydrochloride (3a925)	10 000	mg

ITEMS	AMOUNT	UNIT
Natural or corresponding synthetic chemically define	ed flavourings	
Glycine	8 000	mg
L-Glutamic Acid	5 000	mg
Amino acids, their salts and analogues		
L-Cysteine	6 000	mg
Liquid extract of herbs		
Boldo Plant Extracts	5 000	mg
Artichoke Plant Extracts	5 000	mg
Trace elements		
Iron (E 1)	275	mg
Copper (E 4)	170	mg
Zinc (E 6)	90	mg
Other zootechnical additives		
Dimethylglycine sodium salt	7 000	mg
Carrier		
Glucose syrup, Glycerol, Fructose syrup, Propylene Gly	ycol,	
Xanthan Gum, Maltodextrin, Demineralized water u	ip to Qsp	1 Kg

INDICATION: easily assimilated booster of energy necessary for efforts and protector of the muscles during intense exercise

TARGET ANIMALS: Horses, Camels

USE AND DOSAGE:

Horses: Feed 1 syringe/day for 3-5 days. Begin feeding 3 days prior to event. On the day of the race, feed approximately five to six hours prior to event

Camels: Feed 1 syringe before last training then ½ syringe daily. On the day of the race, feed ½ syringe approximately five to six hours prior to event

STORAGE: Store product in dry place at temperature less than 30 °C, away from sun light and out of the reach of children.

SHELF LIFE: Two years

PACKING: Syringe 50 ml (Net weight 60 g) Country Of Origin: ITALY - EUROPEAN UNION





